

Freigegebene Reinigungs und Desinfektionsmittel

Legende:

1 = Reiniger

2 = Desinfektionsmittel

3 = geeignet für Ultraschallbäder

Hersteller/ Manufacturer/ Fabricante	Handelsname/ Commerical/ Nombre comercial	Typ/ Type/ Tipo
Diverse	NaOH*(1 mol/l, max. 1h!)	2
Acto GmbH	Actosed Endo Terra	2
Advances Sterilization Products	Cidezyme	1,3
	Cidex	2
	Cidex OPA	2
	Enzol	1,3
	NU Cidex	2
Alkapharm	Peralkan	2
Alpro Dental-Produkte GmbH	BIB forte	1,2,3
Anios	Aniosyme P.L.A.	1,2,3
	Anioxyde 1000	2
	Octanius basique	1,2
Antiseptica chem.pharm.Prod. GmbH	Triacid N	1,2,3
Bandelin electronic GmbH	Stammopur DR	1,2,3
	Stammopur DR 8	1,2,3
	Stammopur GR	1,3
	Stammopur R	1,3
B.Braun Medical AG	Stabimed	2
Bode Chemie GmbH	Aseptisol	2,3
	Korsolex AF	2,3
	Korsolex Basic	2,3
Borer Chemie AG	Deconex 36 Intensiv	1,3
	Deconex 50 FF	2
	Deconex 53 Plus	2
	Superborat	1
Deppe, Laboratorium Dr. rer. nat.	EndoStar	2
	InstruPlus	1,2,3
	InstruPlus N	1,2,3
	InstruStar	1,2,3
Ecolab GmbH & Co. OHG	InstruZym	2,3
	Sekusept easy	2
	Sekusept extra N	1,2,3
	Sekusept forte	2,3
	Sekusept Plus	1,2,3
	Sekusept Pulver Classic	1,2,3
Esteer GmbH	Sekusept Pulver Classic + Aktivator	1,2,3
	Ultrademit Konz.	2,3
Fresenius AG	Afid	1,2
	Afid plus (neu)	1,2
	Sporcid FF	1,2
Indeba Industria E Commercio LTDA	Letahdeido	2
Jose Collado S.A.	Darodor 4000 Liquido	1
	Darodor 9000	2
	Darodor Sinaldehhyd 2000	2
Lysoform Dr. Rosemann GmbH	Aldasan 2000	1,2
	Almyrol	1,2
	Desoform	1,2

	Lysiformin 3000	1,2
Merz + Co.	Edisonite Super	3
	Mucadont-IS	2,3
	Mucadont-Zymaktiv	1,3
	Mucocid-T	1,2,3
Orochemie	A 10	2
	ID 210	2
Promeagent AB	Wavacide	2
Schülke & Mayr GmbH	Gigasept	2,3
	Gigasept FF	2,3
	Lyssetol AF	1,2,3
	Lyssetol V	2
Schumacher, Dr. GmbH	Descoton forte	2
	Perfektan Endo	2
Steris	EnzyCare 2	1,3
Weigert, Dr. GmbH & Co.	neodisher LM 2	1
	neodisher mediclean	1
	neodisher mediclean forte	1,3
	neodisher medizym	1,3
Whiteley Industries PTY. Ltd	Aidal Plus	2