

Freigegebene Reinigungs und Desinfektionsmittel

Legende:

1 = Reiniger

2 = Desinfektionsmittel

3 = geeignet für Ultraschallbäder

Hersteller/ Manufacturer/ Fabricante	Handelsname/ Commerical/ Nombre comercial	Typ/ Type/ Tipo
BHT Hygiene Technik GmbH	BHT Scope Cleaner	1
	BHT Scope Desinfectant	2
Bode Chemie GmbH	Dismoclean 24 Vario	1
	Dismoclean 28 alka one	1
	Korsolex-Endo-Cleaner	1
	Korsolex-Endo-Disinfectant	2
Bohrer Chemie AG	Deconex 23 Neutrazym	1
	Deconex 28 Alka One	1
	Deconex 55 endo	2
Deppe, Laboratorium Dr. rer. nat.	Endomat Plus	1,2
DiverseyLever	Sumotox E	1
Ecolab GmbH & Co. OHG	Sekumatic FD	2
	Sekumatic FRE	1
	Sekumatic FR (only in combination with Sekumatic FNZ)	1
	Sekumatic FNZ	4
IMS	Adaspor Ready to Use	2
Medisafe UK, Ltd.	3E-Zyme/HS-Zyme	1
Merz+Co.	Mucapur ED	2
	Mucapur ER	1
Ruhof Corporation	Endozime AW	1
Schülke & Mayr GmbH	Thermosept ED	2
	Thermosept ER	1
	Thermosept RKN-zym	1
	Thermosept DK	2
Schumacher, Dr. GmbH	Thermoton Endo	2
Steris	Instru-Klenz	1,3
Weigert, Dr. GmbH & Co.	neodisher FA	1
	neodisher Fa forte	1,3
	neodisher mediclean	1
	neodisher mediclean forte	1
	neodisher mediklar	5
	neodisher medizym	1
	neodisher Septo DA	2
neodisher Septo DN	2	